Purpose

Structural element used to make endoskeletal prostheses of the lower limbs


- STARFLEX range : Weight limit : 100 kg
- STARMAX range : Tested in accordance with ISO 10328 P6 136 kg
- TITANIUM range : Weight limit : 136 kg

Adjustments

Tilt amplitude : +/- 9° compared to the reference axis.

Set the desired tilt using the 4 screws, in accordance with the limits set.

Spares parts

		BA320STFX	BA320STMX	BA320TI
1	Flat point screw M8	V03I-0816	V03Z-0816	V03I-0812

Restriction on use / Precautions for use

Tightening torques (use key (Ref. : ZOC102))

Screw	Torques (Nm)	Screw	Torques (Nm)
V03I-0816 and V03I-0812	13	V03Z-0816	15

- Assemble the screws using normal (Ref. : Z31-243) locking compound
- The product must be used within the recommended weight limits and in accordance with the assessment made by an approved ortho-prosthetist depending on the patient concerned
- Product maintenance is subject to an assessment made by an approved ortho-prosthetist depending on the conditions of use
- All assembly and adjustment work must be carried out by qualified personnel

Cleaning / Storage

- Stock the product in a dry place, in its original packing (on which the batch number is marked)
- The frequency of inspections to check the equipment is subject to the assessment made by the ortho-prosthetist depending on the patient concerned
- Use a damp cloth to clean the product

Guarantee

The guarantee is valid for one year as from the date of purchase of the product, if the product has been used in accordance with the manufacturer's instructions.

N.B. In the event of a complaint, please enclose a return form correctly filled in, incomplete forms cannot be dealt with.

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