Fitting

- ☞ Grease these contacts area (Réf. : Z10)


To optimize your assembly work, we recommend using these parts in conjunction with other parts of the same class in the JANTON range.

Adjustments

When making the orthosis, it is highly advisable to optimize joint alignment to reduce wear of the contact surfaces.

NOTE : Control joint parallelism in order to reduce wear of surfaces in contact.

Spares parts

		L23	L24
1	Spring	L2307	L2407
2	Piston	L2304	L2304
3	Cap	L2305	L2305
4	Tip	L2306	L2306
5	Axis	L2308	L2308
6	Spindle	L2302	L2402
	Spindle pipe	YA23-0811 Ø 8 x 11 mm	YA23-0609 Ø 6 x 9 mm

Restriction on use / Precautions for use

- ☞ Do not heat the joints
- ☞ The product must be used within the recommended weight limits and in accordance with the assessment made by an approved ortho-prosthetist depending on the patient concerned
- ☞ Product maintenance is subject to an assessment made by an approved ortho-prosthetist depending on the conditions of use
- ☞ All assembly and adjustment work must be carried out by qualified personnel

Cleaning / Storage

- ☞ Stock the product in a dry place, in its original packing (on which the batch number is marked)
- ☞ The frequency of inspections to check the equipment is subject to the assessment made by the ortho-prosthetist depending on the patient concerned

Guarantee

The guarantee is valid for one year as from the date of purchase of the product, if the product has been used in accordance with the manufacturer's instructions.

N.B. In the event of a complaint, please enclose a return form correctly filled in, incomplete forms cannot be dealt with.

EC markings

Groupe G2M hereby declares, in its quality as the manufacturer, that all its parts comply with directive 93/42/EC.

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